

# D-LITE

Specialty scalers and curettes

***D-Lite scalers & curettes are specially designed with you in mind. Achieve peak performance and maximize patient care with durable, color-coded, 100% stainless steel instruments.***

## INNOVATIVE TIP DESIGNS



- **Ion charged technology** increases edge retention
- **Computerized hydraulic forming** ensures 100% angle accuracy.
- Made in America

## D-LITE HANDLE



- #6 Light easy grip handle
- Smooth transition from shank to handle
- Available single & double ended
- 440A stainless steel
- Diamond knurl for optimal grip


## TAKE ACTION, organize your operatory

### Chart your course to patient success


*with strategic and effective instrument maintenance solutions:*

- Color chart your instruments and brighten your tray
- Color coded handles correspond with Gracey & Scaler/Curette mouth charts below.

#### Gracey's



#### Scalers & Curettes



Follow us on Facebook and LinkedIn

1-877-284-8261

www.atitan.com