



DISPOSABLE BONE SCRAPERS & COLLECTORS

BG200

Recommended for ramus, symphysis, and tuberosity areas of the mouth.



BG205

Recommended for the external oblique ridge. Hold at a 35 degree angle on the ridge.



- Sterile Disposable bone scraper
- Always has a fresh blade
- 360* blade is marked in 4 sections
- Expose new cutting surface with the twist of any .050 hex driver

BC100

Osseous Collector, pack of 10



- Sterile Disposable bone collector
- Harvest valuable autogenous bone
- Correct small defects around implant sites and periodontal comprised teeth
- Shut off valve cuts off suction eliminating the chance of drying out bone
- Includes 2 filters
- Sold in packs of 10
- Single patient use only

OTA PERIOTOMES



OTA400 Set of 4
Includes one of each:
P301 Blue Periotome, P302 Gray Periotome, P303 Gold Periotome and P304 Green Periotome.

- Ideal for atraumatic extractions in the molar region
- Available in 2 blade configurations
- Color coded for easy identification

P405



P406



- OTA periotomes available in four different blade configurations - wide, narrow and two angled
- Color coded for easy identification

For redemption of goods in association with a current promotion, please fax or email a copy of your invoice to 716-667-9212 or orders@atitan.com. Please indicate which item codes are to be sent per promotion.